

# X-GEL-70

## EXTRA POWER INSTANT HAND SANITIZER

Germs remain active on hard surfaces for hours, and sometimes even days. These germs are passed on from one person to another when we handle everyday things like telephones, computer keyboards and mice, door handles, hand rails, money, etc. Experts acknowledge that approximately 80% of all common infections are spread by the hands this way.

Opus X-GEL-70, which is formulated using 70% ethanol, is an extra effective weapon in the ever more serious battle against pathogenic bacteria.

Opus X-GEL-70 has been developed to meet the requirements of the care sector for a completely effective, yet quick-acting and nontainting hand steriliser. Combining 70% alcohol gives proven bactericidal effectiveness.

**Added glycerine leaves hands soft, supple, non-sticky and moisturised - does not dry out the skin and leaves no after-odour.**

These advanced properties make Opus X-GEL-70 equally suitable for many medical and care service applications - in fact, any application where there is a risk of cross-infection from hand contact.

Independent laboratory tests prove 99.9999% kill of gram positive and gram negative bacteria.

Complies with EU Approval EN1276 for food use and is effective against *Pseudomonas Aeruginosa*, *Escherichia Coli* ATCC 10536, *Staphylococcus Aureus* ATCC 6538, and *Enterococcus hirae* after a contact time of 5 minutes.

Complies with EU Approval EN13727 for effective use against MRSA, *Salmonella Typhimurium*, *Listeria Monocytogenes* in medical areas.

Apply a small quantity to clean, dry hands and rub well in, paying particular attention to fingernails and creases in the skin. Within 15 to 20 seconds, X-GEL-70 will have evaporated, leaving the skin soft and dry. It may also be used in the same way to sterilise surgical gloves.

**MSDS Data:** An EC safety data sheet providing additional information is available on request.

**Packaging Sizes:**

- 24 x 80ml pump or flip-top
- 12 x 250ml pump or flip-top
- 10 x 500ml pump bottles
- 6 x 1 litre pump bottles
- 6 x 1 litre cartridges for use in dedicated wall-mounted dispenser
- 4 x 5 litre bulk-fill jerrycans

**X-GEL-70 IS USED IN:**

- Clinics
- Dental surgeries
- Doctors' surgeries
- Food preparation areas
- Hospitals
- Kindergartens
- Nursing homes
- Schools
- Colleges
- Butchers' shops
- Veterinary surgeries
- Local Councils
- Supermarkets
- Bakeries



Excalibur Hygiene Limited - Opus Handcare Division  
 Unit 42 Enterprise Trading Estate | Pedmore Road  
 Brierley Hill | West Midlands | DY5 1TX | United Kingdom  
 Tel: +44 (0)1384 671505 | Fax: +44 (0)1384 671256  
 Email: [help@excalibur-hygiene.co.uk](mailto:help@excalibur-hygiene.co.uk)  
 Website: [www.hand-cleaners.co.uk](http://www.hand-cleaners.co.uk)

